


# Removing samples from the Bruker 500 MHz autosampler



McGILL CHEMISTRY  
CHARACTERIZATION


- Check which samples can be removed

*To bring up the sample positions display*



*Green and yellow positions can be removed*

*From the sample changer to the same position number in any of the three racks*



*Older samples: from a numbered position to no-number positions in the front*



- Robin will remove samples from the no-number positions and put them in her office for at least a week, then will give them away
- *If you're looking for your sample, it should either be in the position number you left it in, the same position number in any of the three racks, the front positions where there are no numbers, or in Robin's office – ask her!*
- *You can take samples from the "FREE" rack, but do not put any samples there*